

LEITZ SM-LUX MICROSCOPE



DESCRIPTION

We offer the Leitz SM-LUX microscope for sale, which is in excellent physical and working condition. This optical instrument by Ernst Leitz Wetzlar has been recently overhauled to its original factory specifications. Featuring a WF 10x microscope eyepiece and quintuple nosepiece, with coaxial coarse and fine adjustment focusing, and a built-in light source, the Leitz Wetzlar SM-LUX microscope makes viewing samples easy.

The microscope objectives of 10x, 40x, and 100x are all Leitz-quality optics. A Flip Top Abbe Condenser 1.25 N.A. comes standard on this Leitz SM-LUX microscope model.

PRODUCT SPECIFICATIONS

Style	Binocular
Eyepiece	WF 10x
Nosepiece	Quintuple
Focusing	Coaxial coarse & fine adjustment
Stage	Mechanical stage
Objectives	Leitz optics 10x, 40x, 100x
Sub-stage	Flip Top Abbe Condenser 1.25 N.A.
Light Source	Built in light source with variable light intensity
Warranty	One Year